

ICU X-Ray Protective Curtains



detail text for image:
ICU X-Ray Protective Curtains



AmRay Protection

- Bespoke design solutions
- Standard curtain height of 200cm
- Heavy duty track to support curtain weight and provide ease of closing and opening and storage
- Ceiling mounted curtain track (Model XC01)
- Fold-away curtain solution allows curtain to be stored away against the wall (Model FXC01)
- Robust design with performance eyelets
- PU coated material with lead inner core
- Colour options

AmRay Protection

- Bespoke design solutions
- Standard curtain of 200cm
- Heavy duty track to support curtain weight and provide ease of closing, opening and storage
- Ceiling mounted and fold-away curtain solution allows curtain to be stored away against the wall (Model XC01)
- Robust design with performance eyelets
- PU coated material with lead inner core



- Colour options

AmRay's **ICU X-ray Protective Curtains** provide protection to staff, patients and the general public. Can be conveniently tied-back when not in use.