




X-Ray Protection Control Booth



detail text for image:
X-Ray Protection Control Booth



AmRay Protection

- Bespoke Design Service to meet needs according to use

AmRay Protection

- Safe and optimal viewing environment for X-ray procedures
- Bespoke design service to meet X-ray viewing needs
- Work station housed within cabin if required
- Lead core panels faced on both sides
- Lead glass viewing panels incorporated within the design
- Aluminium frame
- Choice of colour finish

AmRay's **X-ray Protection Control Booths** provide the user with a safe viewing environment during X-ray procedures. The cabins are constructed from panels with a lead core that are faced on both sides. Lead glass viewing panels are incorporated into a design with the complete structure being housed within an aluminium frame. Once assembled the whole system becomes a self-supporting unit.

